


ANKLE

LOCKING PLATE & SCREW SYSTEM

SYSTEM OVERVIEW



ANKLE LOCKING PLATE & SCREW SYSTEM SYSTEM OVERVIEW

The advanced principle of TDM USA is to positively impact the healthcare system by partnering with medical professionals to define and deliver value-generating solutions. The Ankle Locking Plate and Screw System is a component of the extremity fixation product family that constitutes the **ADVANCED PRINCIPLE IMPLANT SYSTEMS**.



THE ANKLE LOCKING PLATE AND SCREW SYSTEM offers multiple fixation options in a variety of configurations. Plates in this system are intended for use in the distal fibula, distal tibia, and calcaneus. The low profile plates are constructed of titanium and are used with solid locking and non-locking screws. The titanium alloy screws range in length from 10-50mm and are available in 3.5mm diameter for proximal placement and 2.7mm diameter for distal placement.

ANKLE LOCKING PLATE & SCREW SYSTEM

IMPLANT ORDERING INFORMATION

FIBULAR LOCKING PLATE, LATERAL

Fibular Locking Plate, Lateral, Right, 3H	61mm	166-13203
Fibular Locking Plate, Lateral, Right, 4H	72mm	166-13204
Fibular Locking Plate, Lateral, Right, 5H	83mm	166-13205
Fibular Locking Plate, Lateral, Right, 7H	105mm	166-13207
Fibular Locking Plate, Lateral, Left, 3H	61mm	166-13103
Fibular Locking Plate, Lateral, Left, 4H	72mm	166-13104
Fibular Locking Plate, Lateral, Left, 5H	83mm	166-13105
Fibular Locking Plate, Lateral, Left, 7H	105mm	166-13107



FIBULAR LOCKING PLATE, POSTERIOR

Fibular Locking Plate, Posterior, Right, 3H	67mm	166-14203
Fibular Locking Plate, Posterior, Right, 4H	78mm	166-14204
Fibular Locking Plate, Posterior, Right, 5H	88mm	166-14205
Fibular Locking Plate, Posterior, Left, 3H	67mm	166-14103
Fibular Locking Plate, Posterior, Left, 4H	78mm	166-14104
Fibular Locking Plate, Posterior, Left, 5H	88mm	166-14105



DISTAL TIBIA LOCKING PLATE, ANTERIOR

Distal Tibia Locking Plate, Anterior, Right, 3H	55mm	160-11203
Distal Tibia Locking Plate, Anterior, Right, 4H	62mm	160-11204
Distal Tibia Locking Plate, Anterior, Right, 5H	71mm	160-11205
Distal Tibia Locking Plate, Anterior, Right, 6H	79mm	160-11206
Distal Tibia Locking Plate, Anterior, Left, 3H	55mm	160-11103
Distal Tibia Locking Plate, Anterior, Left, 4H	62mm	160-11104
Distal Tibia Locking Plate, Anterior, Left, 5H	71mm	160-11105
Distal Tibia Locking Plate, Anterior, Left, 6H	79mm	160-11106



DISTAL TIBIA LOCKING PLATE, POSTERIOR

Distal Tibia Locking Plate, Posterior, Right, 3H	60mm	161-14203
Distal Tibia Locking Plate, Posterior, Right, 4H	71mm	161-14204
Distal Tibia Locking Plate, Posterior, Right, 5H	82mm	161-14205
Distal Tibia Locking Plate, Posterior, Right, 6H	93mm	161-14206
Distal Tibia Locking Plate, Posterior, Left, 3H	60mm	161-14103
Distal Tibia Locking Plate, Posterior, Left, 4H	71mm	161-14104
Distal Tibia Locking Plate, Posterior, Left, 5H	82mm	161-14105
Distal Tibia Locking Plate, Posterior, Left, 6H	93mm	161-14106



DISTAL TIBIA LOCKING PLATES, MEDIAL

Distal Tibia Locking Plate, Medial, 3H	58mm	161-15003
Distal Tibia Locking Plate, Medial, 4H	68mm	161-15004
Distal Tibia Locking Plate, Medial, Hook, 3H	60mm	161-16003
Distal Tibia Locking Plate, Medial, Hook, 4H	70mm	161-16004



3.5mm LOCKING CALCANEUS PLATE

3.5mm Locking Calcaneus Plate	50mm	144-10050
3.5mm Locking Calcaneus Plate	60mm	144-10060
3.5mm Locking Calcaneus Plate	70mm	144-10070



ANKLE LOCKING PLATE & SCREW SYSTEM

IMPLANT ORDERING INFORMATION/INSTRUMENT OVERVIEW

2.7mm LOCKING SCREW

10-40mm

227-260xx*

2.7mm CORTICAL SCREW, LOW PROFILE HEAD

10-40mm

227-270xx*

3.5mm LOCKING SCREW

10-50mm

235-260xx*

3.5mm CORTICAL SCREW

10-40mm

235-210xx*

(All screws in 2mm increments)

*xx = Screw length



Temporary Fixation K-wire 1.4mm 901-16014 	Plate Bender, Small Left Right 901-06035 901-06036 
Drill Bit 2.0 X 130mm (Yellow) 901-01220 2.5 X 130mm (White) 901-01125 2.7 X 130mm (Blue) 901-22127 3.5 X 160mm (Green) 901-22135 	Drill Bit, AO QC 2.0 X 130mm (Yellow) 901-31120 2.5 X 130mm (White) 901-31025 2.7 X 130mm (Blue) 901-31127 3.5 X 160mm (Green) 901-31135 
Locking Drill Sleeve 2.0mm (Yellow) 901-12520 2.7mm (Blue) 901-12527 	Locking Drill Sleeve, Single Handle 2.0mm (Yellow) 901-02220 2.7mm (Blue) 901-02227 
Drill Sleeve, Single Handle 2.5mm (White) 901-13325 3.5mm (Green) 901-13335 	Depth Gauge 2.7/3.5/4.0mm 901-04027 
Screwdriver, Self-Retaining 2.0mm Hex 901-35020 2.5mm Hex 901-35025 	Torque Limiting Handle, AO QC 0.8Nm 901-17008 1.5Nm 901-17015 
Screwdriver Shaft, Self-Retaining, AO QC 2.0mm Hex 901-15420 2.5mm Hex 901-15425 	Forceps 901-08001 
Drill Guide Block, Distal Tibia Plate Anterior Right: 901-22269 Left: 901-22268 Posterior Right: 901-23269 Left: 901-23268	Drill Guide Block, Fibular Plate Lateral Right: 901-20269 Left: 901-20268 Posterior Right: 901-21269 Left: 901-21268

*See plate image on previous page

*See plate image on previous page